

Open Existing Project[Name 1](#) [Name 2](#) [Name 3](#) Log in to have unlimited
access to your projects.☒ Disable cookies [Explain](#)Enter a DNA sequence, or select from other options, to identify cut sites.
Once you submit a sequence, you may choose to customize your digest.**1. Input or choose sequence.** [Explain](#)**Paste Sequence**Plain or FASTA/
Max: 300 KBasesOR **Upload File**
(xx MB max.)

No file selected.

OR **Enter Gen Bank**[Browse Gen Bank](#)OR **Select Sequence****2. Update preferences.** [Explain](#)☒ Circular [Additional Preferences](#)**3. Name project (optional).** [Explain](#)

Additional Preferences

Enzymes to include:

- ☒ NEB enzymes (Type II & Type III)
 - ☐ [Time-Saver™](#) qualified enzymes
 - ☐ All commercially available specificities
 - ☐ All specificities
 - ☐ All & defined oligonucleotide sequences
 - ☐ Only defined oligonucleotide sequences
- [Define Oligos](#)

In addition to Type II & commercially-available Type III enzymes, also include:

- ☒ Type I & III enzymes
- ☒ Homing endonucleases
- ☒ Nicking enzymes

Methylation sensitivities to include:

- ☒ CpG methylation
- ☒ Dam methylation
- ☒ Dcm methylation
- ☒ EcoBI methylation
- ☒ EcoKI methylation

Minimum ORF length:

 aa

Genetic code for ORFs:

⌵ Select Genetic Code

Process this region only: – bp

Save

Cancel

3. Name project (optional). [Explain](#)

Save

Submit

Standard Digest for PURExpr prmr

[Make PDF](#)

[Print](#)

Use the menu to explore this digest, or to create a custom digest.

Graphical View

Enzyme List

Sequence

ORF Summary

Customize Digest

Results For:


- Linear
- NEB single cutter restriction enzymes
- Main non-overlapping, min. 100 aa ORFs

Sequence Information:

2817 bp
GC = 34% AT = 66%

[Update Preferences](#)



[Start Over](#)






Display

- ☒ 1 cutters
- ☒ 2 cutters
- ☒ 3 cutters

Cleavage Code

-  Blunt-End Cut
-  Cuts 1 strand
-  5' Extension
-  3' Extension

Enzyme Code

-  NEB
-  Other supplier
-  Not commercially available

[Update Colors](#)

Methylation

- * CpG methylation
- # dam/dcm methylation
- () Ambiguous site

Standard Digest for PURExpr prmr

Make PDF

Print

Graphical View

Enzyme List

Sequence

ORF Summary

Customize Digest

Results For:

- Linear
- NEB single cutter restriction enzymes
- Main non-overlapping, min. 100 aa ORFs

Sequence Information:

2817 bp
GC = 34% AT = 66%

Update Preferences

Start Over

Cut Position ▼	Enzyme ▲	Recognition Sequence	Me Sens.
2	MluCI	^AATT_	
5	AseI	AT^TA_AT	
5	MseI	T^TA_A	
6	FspEI	CC(N)12^NNNN_	
6	PleI	GAGTCNNNN^N_	
6	MlyI	GAGTC(N)5	
7	FspEI	CC(N)12^NNNN_	
7	MspJI	CNNR(N)9^NNNN_	
12	HinfI	G^ANT_C	CpG
18	SfcI	C^TRYA_G	
19	FspEI	CC(N)12^NNNN_	CpG
19	BsmAI	GTCTCN^NNNN_	
19	BsaI	GGTCTCN^NNNN_	
19	BcoDI	GTCTCN^NNNN_	
26	MspJI	CNNR(N)9^NNNN_	
36	HpyCH4III	AC_N^GT	
42	FspEI	CC(N)12^NNNN_	
42	MspJI	CNNR(N)9^NNNN_	
42	MspJI	CNNR(N)9^NNNN_	
44	XbaI	T^CTAG_A	

Display

- ☒ 0 cutters
- ☒ All sites
- ☒ 2 cutters
- ☒ Flanking enzymes



NEB Restriction Enzyme
Activity/Performance Chart

Standard Digest for LITMUS39

[Make PDF](#)

[Print](#)

Graphical View

Enzyme List

Sequence

ORF Summary

Customize Digest

Results For:

- Circular
- NEB single cutter restriction enzymes
- Main non-overlapping, min. 100 aa ORFs

Sequence Information:

2817 bp
GC = 34% AT = 66%

[Update Preferences](#)

[Start Over](#)

```
1  GTTAACTACG TCAGGTGGCA CTTTTCGGGG AAATGTGCGC GGAACCCCTA TTTGTTTATT
61  TTTCTAAATA CATTCAAATA TGTATCCGCT CATGAGACAA TAACCCTGAT AAATGCTTCA
121 ATAATATTGA AAAAGGAAGA GTATGAGTAT TCAACATTTT CGTGTCGCCC TTATTCCCTT
181 TTTTGCGGCA TTTTGCCTTC CTGTTTTTGC TCACCCAGAA ACGCTGGTGA AAGTAAAAGA
241 TGCTGAAGAT CAGTTGGGTG CACGAGTGGG TTACATCGAA CTGGATCTCA ACAGCGGTAA
301 GATCCTTGAG AGTTTTTCGCC CCGAAGAACG TTCTCCAATG ATGAGCACTT TTAAAGTTCT
361 GCTATGTGGC GCGGTATTAT CCCGTGTTGA CGCCGGGCAA GAGCAACTCG GTCGCCGCAT
421 ACACTATTCT CAGAATGACT TGGTTGAGTA CTCACCAGTC ACAGAAAAGC ATCTTACGGA
481 TGGCATGACA GTAAGAGAAT TATGCAGTGC TGCCATAACC ATGAGTGATA ACACTGCGGC
541 CAACTTACTT CTGACAACGA TCGGAGGACC GAAGGAGCTA ACCGCTTTTT TGCACAACAT
601 GGGGGATCAT GTAACTCGCC TTGATCGTTG GGAACCGGAG CTGAATGAAG CCATACCAAA
661 CGACGAGCGT GACACCACGA TGCCTGTAGC AATGGCAACA ACGTTGCGCA AACTATTAAC
721 TGGCGAACTA CTTACTCTAG CTTCCCGGCA ACAATTAATA GACTGGATGG AGGCGGATAA
781 AGTTGCAGGA CCACTTCTGC GCTCGGCCCT TCCGGCTGGC TGGTTTATTG CTGATAAATC
841 TGGAGCCGGT GAGCGTGGGT CTCGCGGTAT CATTGCAGCA CTGGGGCCAG ATGGTAAGCC
901 CTCCCGTATC GTAGTTATCT ACACGACGGG GAGTCAGGCA ACTATGGATG AACGAAATAG
961 ACAGATCGCT GAGATAGGTG CCTCACTGAT TAAGCATTGG TAACTGTCAG ACCAAGTTTA
1021 CTCATATATA CTTTAGATTG ATTTACCCCG GTTGATAATC AGAAAAGCCC CAAAAACAGG
1081 AAGATTGTAT AAGCAAATAT TTAAATTGTA AACGTTAATA TTTTGTTAAA ATTCGCGTTA
1141 AATTTTTGTT AAATCAGCTC ATTTTTTAAC CAATAGGCCG AAATCGGCAA AATCCCTTAT
1201 AAATCAAAAG AATAGCCCGA GATAGGGTTG AGTGTTGTTC CAGTTTGGA CAAGAGTCCA
1261 CTATTAAAGA ACGTGGACTC CAACGTCAAA GGGCGAAAAA CCGTCTATCA GGGCGATGGC
1321 CCACTACGTG AACCATCACC CAAATCAAGT TTTTGGGGT CGAGGTGCCG TAAAGCACTA
1381 AATCGGAACC CTAAAGGGAG CCCCCGATTT AGAGCTTGAC GGGGAAAGCG AACGTGGCGA
1441 GAAAGGAAGG GAAGAAAGCG AAAGGAGCGG GCGCTAGGGC GCTGGCAAGT GTAGCGGTCA
1501 CGCTGCGCGT AACCACCACA CCCGCCGCGC TTAATGCGCC GCTACAGGGC GCGTAAAAGG
1561 ATCTAGGTGA AGATCCTTTT TGATAATCTC ATGACC AAAA TCCCTTAACG TGAGTTTTTCG
1621 TTCCACTGAG CGTCAGACCC CGTAGAAAAG ATCAAAGGAT CTTCTTGAGA TCCTTTTTTTT
1681 CTGCGCGTAA TCTGCTGCTT GCAAACAAAA AAACCACCGC TACCAGCGGT GGTGTTGTTG
1741 CCGGATCAAG AGCTACCAAC TCTTTTTCCG AAGGTA ACTG GCTTCAGCAG AGCGCAGATA
```



Standard Digest for LITMUS39

Make PDF

Print

Graphical View

Enzyme List

Sequence

ORF Summary

Customize Digest

Results For:

- Circular
- NEB single cutter restriction enzymes
- Main non-overlapping, min. 100 aa ORFs

Sequence Information:

2817 bp
GC = 34% AT = 66%

Update Preferences

Start Over

Gene	Product	a.a. Seq.	Coordinates	Protein ID	Closest Enzyme		Flanking Enzymes
					5´ End	3´ End	
bla	beta-lactamase	286 aa	143..1003	–	EarI	BsaJI	Show
lacZalpha	beta-galactosidase alpha fragment	121 aa	2424..2789	–	*BsrBI	HpaI	Show

Custom Digest for pUC19

[Make PDF](#)

[Print](#)

Graphical View

Enzyme List

Fragments

Gel

Digested With:

- EcoRI
- HindIII

Sequence Information:

2686 bp
GC = 51% AT = 49%

[Update Preferences](#)

[Start Over](#)

EcoRI

[NEB #R0101](#)

Features:

- 5 minute Time-Saver
- Buffer: NEBuffer EcoRI
- Salt: 50 mM NaCl
- Main: 100 mM Tris-HCl
- pH: 7.5
- Mg: 10 mM MgCl₂
- Triton: 0.025%
- Rxn temp.: 37 °C
- Neoschizomers: None
- Isoschizomers: None
- Sites in sequence: [1](#)
- End produced at 396: 5' overhang: AATT

[REBASE page](#)



**NEB Restriction Enzyme
Activity/Performance Chart**

Methylation:

! This site overlaps a CpG site: GAATTCg ,
overlapping sites may impair cleavage if methylated.

Enzymes producing compatible ends:

Enzymes with no sites in pUC19:

Create Custom Digest for pUC19

Only enzymes that cut the sequence and match your selection of availability are shown. Non-NEB neoschizomers are not listed, but you can enter it in the search field and its prototype will be selected.

Jump to enzyme(s) or oligo sequence(s) 

Select	Enzyme ▼	Recognition Sequence	Cuts
All			
<input type="checkbox"/>	AatII	GACGTC	1
<input type="checkbox"/>	Acc65I	GGTACC	1
<input type="checkbox"/>	AccI	GTMKAC	1
<input type="checkbox"/>	AciI	CCGC	34
<input type="checkbox"/>	AclI	AACGTT	2
<input type="checkbox"/>	AclI	CTGAAG(N)14NN	2
<input type="checkbox"/>	AfiIII	ACRYGT	1
<input type="checkbox"/>	AhdI	GACNNNNNGTC 1	1
<input type="checkbox"/>	AluI	AGCT	16
<input type="checkbox"/>	AlwI	GGATCNNNNN	10
<input type="checkbox"/>	AlwNI	CAGNNNCTG	1
<input type="checkbox"/>	ApaLI	GTGCAC	3
<input checked="" type="checkbox"/>	ApeKI	GCWGC	12
<input type="checkbox"/>	ApoI	RAATTY	1
<input type="checkbox"/>	AseI	ATTAAT	3
<input type="checkbox"/>	AvaI	CYCGRG	1
<input type="checkbox"/>	AvaII	GGWCC	2
<input checked="" type="checkbox"/>	BaeGI	GKGCMC	3
<input type="checkbox"/>	BamHI	GGATCC	1
<input type="checkbox"/>	BanI	GGYRCC	4

Filter

- ☒ Enzymes with compatible buffers
- ☒ Enzymes producing blunt ends
- ☒ Enzymes producing 5' overhangs
- ☒ Enzymes producing 3' overhangs
- ☒ Enzymes with a particular site length
- ☒ Enzymes cutting N times
- ☒ Non-palindromic enzymes

Selections:

- ApeKI
- BaeGI
- EcoRI

Digest

Reset

Custom Digest for pUC19

[Make PDF](#) [Print](#)

Graphical View

Enzyme List

Fragments

Gel

Digested With:


- EcoRI
- HindIII

Sequence Information:

2686 bp
GC = 51% AT = 49%

[Update Preferences](#)

[Start Over](#)







Display




- ☒ Circular
☐ Linear

To zoom into cut sites,
use linear display.

Cleavage Code

-  Blunt-End Cut
 Cuts 1 strand
 5' Extension
 3' Extension

Enzyme Code

-  NEB
 Other supplier
 Not commercially available

[Update Colors](#)

Methylation

- * CpG methylation
dam/dcm methylation
() Ambiguous site

ORFs

a: bla
b: lacZalpha

Custom Digest for pUC19

[Make PDF](#)

[Print](#)

Graphical View

Enzyme List

Fragments

Gel

Digested With:

- EcoRI
- HindIII


Sequence Information:

2686 bp

GC = 51% AT = 49%

[Update Preferences](#)

[Start Over](#)



Display

- ☐ Circular
- ☒ Linear

To zoom into cut sites,
use linear display.

Cleavage Code

- ✂ Blunt-End Cut
- ▼ Cuts 1 strand
- ▼ 5' Extension
- ▲ 3' Extension

Enzyme Code

- NEB
- Other supplier
- Not commercially available

[Update Colors](#)

Methylation

- * CpG methylation
- # dam/dcm methylation
- () Ambiguous site

ORFs

a: bla

b: lacZalpha

Custom Digest for pUC19

[Make PDF](#)

[Print](#)

Graphical View

Enzyme List

Fragments

Gel

Sort By

#	Enzyme ▼	Recognition Sequence	Cuts	Sites & Flanks	Cut Positions	Blunt/ Overhang
1	EcoRI	G [^] AATT [^] C	1	Show	*396/400	5´
2	HindIII	A [^] AGCT [^] T	1	Show	447/451	5´

[Save as Text](#)

Digested With:

- EcoRI
- HindIII

Sequence Information:

2686 bp
GC = 51% AT = 49%

[Update Preferences](#)

[Start Over](#)

Display

- ☒ 0 cutters
- ☒ All sites
- ☒ 2 cutters
- ☒ Flanking enzymes



**NEB Restriction Enzyme
Activity/Performance Chart**

Custom Digest for pUC19

Make PDF

Print

Graphical View

Enzyme List

Fragments

Gel

Digested With:

- EcoRI
- HindIII

Sequence Information:

2686 bp
GC = 51% AT = 49%

Update Preferences

Start Over

Sort by

#	Ends	Coordinates	Length (bp)
1	HindIII-EcoRI	448-396	26351
2	EcoRI-HindIII	397-447	51

Save as Text


#	Ends	Coordinates	Length (bp)
1	EcoRI-HindIII	397-447	51
2	HindIII-EcoRI	448-396	2635

Custom Digest for pUC19

[Make PDF](#)

[Print](#)

This virtual gel was generated by interpolating experimental data. [See details.](#)



#	Ends	Coordinates	Length (bp)
1	HindIII-EcoRI	448-396	26351
2	EcoRI-HindIII	397-447	51

Graphical View

Enzyme List

Fragments

Gel

Digested With:

- EcoRI
- HindIII

Sequence Information:

2686 bp
GC = 51% AT = 49%

[Update Preferences](#)

[Start Over](#)

[Download Gel](#)

Gel Length

– mm

[Update](#)

Color

- ☒ White on black
- ☐ Black on white

◆ Gel Type

◆ Marker

◆ DNA Type